

## Specular Microscope

## SPM-700

## Specifications

Capturing of corneal endothelial cell	Capturing position	Capturing range	0.25mm×0.55mm (W×H)
		Center	1 point
		Paracenter	6 points(2,4,6,8,10 and 12 o'clock directions)
		Periphery (optic angle: 27 degrees)	10 points (1,2,4,5,6,7,8,10,11 and 12 o'clock directions)
Measurement of corneal thickness	Range of corneal thickness measurement	400 to 750 $\mu$ m (step:1 $\mu$ m)	
Analysis parameter	[Number]	[cells]	Number of endothelial cells
	[CD]	[cell/mm <sup>2</sup> ]	Density of endothelial cells
	[AVG]	[ $\mu$ m <sup>2</sup> ]	Average endothelial cell area
	[SD]	[ $\mu$ m <sup>2</sup> ]	Standard deviation of cell area
	[CV]	[%]	Coefficient of variation of cell area
	[Max]	[ $\mu$ m <sup>2</sup> ]	Max.cell area
	[Min]	[ $\mu$ m <sup>2</sup> ]	Min.cell area
Histogram	[6A]	[%]	Rate of cell hexagonality
	Polymegathism		
Pleomorphism			
Monitor	10.4 inch touch panel colored LCD(XGA)		
Printer	Thermail printer (paper width 58mm)		
External interface	USB-A × 2, USB-B × 1, Ethernet (10/100 Mbps) × 1		
Source voltage / frequency	AC100V-240V,50/60 Hz		
Power consumption	90VA		
Power saving function	OFF,3,5,10min(switchable)		
Size	H(503mm) × W(271mm) × D(459mm)		
Weight	19kg		

SPM-700  
Standard Accessories

- Operation manual
- Power cord
- Printer paper
- Fuse
- Dust cover
- Chinrest paper
- Chinrest paper pin

Design and specifications are subject to change without notice.

Manufacturer \_\_\_\_\_

**Rexxam**

Rexxam Co.,Ltd.

958 Ikeuchi, Konan, Takamatsu-shi,  
Kagawa 761-1494 Japan

Contact \_\_\_\_\_

MEC Sales Division  
2-8-4 Kandatsukasa-machi, Chiyoda-ku, Tokyo 101-0048 Japan  
TEL +81-3-3256-7701 FAX +81-3-3256-7702  
E-mail: eye@rexam.co.jp URL: http://www.rexxam.co.jp

Distributed by \_\_\_\_\_



Specular Microscope

Specular microscope with easy  
operation and speedy analysis

**SPM-700**

# Specular microscope with easy operation and speedy analysis

## Easy, Speedy and Accurate

### Quick Measurement & Analysis

By simple touch panel operation, alignment is achieved automatically, images are captured continuously in 2 sec., and analyzed in 1 sec.  
High speed and accuracy specular microscope has been realized.



#### Multiple Measurement Points

Total 17 measurement points including center, 6 in the paracenter and 10 in the periphery can be measured in the range of 0.25mm by 0.55mm.



#### Edit Function

This function enables to edit the contrast, brightness and analysis result of the endothelial cell image captured. Also, it allows to remove cells, add/delete lines and divide/merge cells.




#### Continuous Capturing of 16 Images

16 images are captured in 2 sec. with our unique zoom function and auto alignment function by touching paracenter area.  
※ Full-auto, semi-auto and manual can be selected in the operation mode.



#### Speedy Analysis Function

After the measurements, the best image is selected automatically from 16 images. After selecting the best image, analysis is finished in 1 sec.  
※ The image can be selected from 16 images manually.



#### Corneal Thickness Measurement

It is possible to capture the endothelial cell and to take a measurement of corneal thickness at the same time.

#### 2 Manual Analyses

There are 2 manual analyses, center method and frame method.



center method



frame method

### Simple & Easy Operation


The monitor can swivel 90 degrees each from center horizontally and tilt 40 degrees upward.  
The swivel/tilt function allows both operator and patient's easy measurement and satisfaction.  
The high-intensity colored LCD with touch panel is equipped.

#### Wide Screen


10.4 inch wide color screen. The swivel/tilt function allows the operator to support easily the patient during operation.



#### Output of Measurement and Analysis Result



Built-in printer output



External report output



#### 4 Types of Display Mode

The display mode can be selected from ① image of endothelial cell, ② trace display, ③ area display and ④ pleomorphism display.



① image of endothelial cell

② trace display



③ area display

④ pleomorphism display

#### Electric Chinrest

It is easy to align the eye position of the patient with the eye mark.

