



# IQ 532™ Laser System

532 nm Laser with DualSense™

## The Advantages of Innovation

- High power with 2500 mW of deliverable laser power
- High speed with pulse durations from 10 - 3000 ms
- DualSense™ provides quick and simple selection of multiple delivery devices for RFID and SMA
- Intuitive graphic touch screen interface with high contrast color display
- Voice confirmation to aid surgical technique
- **Advancing Technology with Optional MicroPulse™ Module**
  - Tissue sparing procedures to delivery energy below levels of standard photocoagulation
  - Reduces thermal spread by confining heat to the target area<sup>1</sup>

## Ergonomic and Easy to Use

- Convenient 3-knob control console
- Programmable memory presets

## Optional Accessories

- **Full-Featured Remote Control**
  - Compact design for easy placement or use in sterile field
  - View displays and adjust parameters from 2 vantage points for increased convenience and efficiency
- **Wireless Footswitch**
  - No cord, no clutter, no limitations
  - Allows precise physician control over power settings
  - Available with power-adjust to control laser actuation and power settings



**IRIDEX**

Elegantly simple solutions™

# IQ 532™ Laser System

## Specifications

<b>Wavelength:</b>	532 nm Green
<b>Weight:</b>	19.2 lb (9.0 kg)
<b>Dimensions:</b>	8.5" H x 12" W x 14" D (21 cm x 30 cm x 35 cm)
<b>Connector Type:</b>	RFID   Resistor
<b>Electrical:</b>	100–240 VAC , 50/60 Hz
<b>Cooling:</b>	Whisper fan
<b>Exposure Duration:</b>	CW-Pulse™: 10–3000 ms
<b>Exposure Interval:</b>	CW-Pulse: 10–3000 ms
<b>MicroPulse™ Duration:</b>	MicroPulse: 0.05–1.00 ms  MicroPulse
<b>MicroPulse Interval:</b>	MicroPulse: 1.00–10.00 ms
<b>Aiming Laser:</b>	Diode laser, 635 nm nominal
<b>Delivery Device</b>	Portable SLA: 0–1800 mW
<b>Power Output:</b>	LIO: 0–2000 mW EndoProbe®: 0–2000 mW

1. Luttrull JK, Sramek C, Palanker D, Spink CJ, Musch DC. Long-term safety, high-resolution imaging, and tissue temperature modeling of subvisible diode micropulse photocoagulation for retinovascular macular edema. *Retina* 2011.

 0086

Specifications are subject to change without notice. EndoProbe, IRIDEX and the IRIDEX logo are registered trademarks and IQ 532, DualSense, MicroPulse and CW-Pulse are trademarks of IRIDEX Corporation. All other trademarks are the property of their respective owners.

Products are covered by one or more of the following U.S. patents: 5,372,595; 5,511,085; 5,521,932; 5,663,979; 5,982,789; 5,979,554; 6,141,143; 6,144,484; 6,222,869; 6,327,291; 6,377,599; 6,540,391; 6,733,490; 6,800,076; 7,537,593; 7,766,904; 7,771,417; and 7,909,816. Other U.S. and international patents pending.

 Emergo Europe  
Molenstraat 15, 2513 BH The Hague, The Netherlands, Tel: (31) (0) 70 345-8570, Fax: (31) (0) 70 346-7299



# IRIDEX

Elegantly simple solutions™



IRIDEX | 1212 Terra Bella Avenue | Mountain View, CA 94043 | 800.388.4747 (U.S. inquiries) | info@iridex.com (U.S. & int'l inquiries) | www.iridex.com

©2011 IRIDEX. Printed in U.S.A. LT0488 15/10/2011