

VT-1000/VT-1380



CE FDA

Specification

Table top size	137x65cm (VT-1000)
	95x52cm (VT-1380)
Table Elevator Height	63-93cm
Maximum load	80kgs
Input	100V-220V
Frequency	50HZ-60HZ
Weight	61KG

Shanghai VisuScience Meditech Co.,Ltd.

Add: No. 344 Sanlin Road, Pudong New Area, Shanghai, China

Tel:+86-21- 34973659 | Fax: +86-21- 34973659

E-mail:info@visuscience.com

Follow us on



www.visuscience.com