

## Ocular Landers ROP Lens Attachment



### Product Code

OI-LROP



*Designed With:  
Maurice B. Landers, III, M.D.  
Chapel Hill, NC*

### Design

- Engraving on rim enables to precisely measure the number of clock hours of Neo-Vascularization when viewing retinopathy of prematurity (ROP)
- Monofilament line provides separation of quadrants for quick orientation.
- Rim of attachment is conveniently marked at hour and half hour locations.
- Fits onto anterior side of Ocular 28D Indirect lenses\*

\*lens design with diamond knurl pattern only

### Cleaning & Disinfecting

Per Cleaning Method 2



2255 116th Ave NE, Bellevue, Washington 98004-3039 USA  
T: 425-455-5200 or 800-888-6616 F: 425-462-6669  
E: [ocular@ocular-instruments.com](mailto:ocular@ocular-instruments.com) I: [www.ocular-instruments.com](http://www.ocular-instruments.com)

© 2006 OCULAR INSTRUMENTS  
10069-3000